

## Product datasheet for **RR206404L3V**

### **Cnga4 (NM\_053496) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Cnga4 (NM_053496) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cnga4
Synonyms:	Cgn2; rOCNC2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_053496
ORF Size:	1725 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206404).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_053496.2</a> , <a href="#">NP_445948.2</a>
RefSeq Size:	2355 bp
RefSeq ORF:	1728 bp
Locus ID:	85258
UniProt ID:	<a href="#">Q64359</a>
Cytogenetics:	1q32
Gene Summary:	olfactory cyclic nucleotide gated cation channel [RGD, Feb 2006]



[View online »](#)