

## Product datasheet for **RR206410L4V**

### Gsta3 (NM\_001009920) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Gsta3 (NM_001009920) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gsta3
Synonyms:	Gsta5; Yc2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001009920
ORF Size:	663 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206410).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001009920.2</a> , <a href="#">NP_001009920.1</a>
RefSeq Size:	1508 bp
RefSeq ORF:	666 bp
Locus ID:	494500
UniProt ID:	<a href="#">P46418</a>
Cytogenetics:	9q13
Gene Summary:	Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles. Has substantial activity toward aflatoxin B1-8,9-epoxide. [UniProtKB/Swiss-Prot Function]



[View online »](#)