

Product datasheet for **RR206485L3V**

Galk2 (NM_001013919) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Galk2 (NM_001013919) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Galk2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001013919
ORF Size:	1374 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206485).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001013919.1 , NP_001013941.1
RefSeq Size:	1797 bp
RefSeq ORF:	1377 bp
Locus ID:	296117
UniProt ID:	Q5XIG6
Cytogenetics:	3q36
Gene Summary:	Acts on GalNAc. Also acts as a galactokinase when galactose is present at high concentrations (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)