

Product datasheet for **RR206493L3V**

Ehd3 (NM_138890) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Ehd3 (NM_138890) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ehd3
Synonyms:	Ehd2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_138890
ORF Size:	1605 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206493).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138890.2 , NP_620245.2
RefSeq Size:	2108 bp
RefSeq ORF:	1608 bp
Locus ID:	192249
UniProt ID:	Q8R491
Cytogenetics:	6q13
Gene Summary:	mouse homolog is involved in regulating vesicular structure microtubule-dependent movement [RGD, Feb 2006]



[View online »](#)