

Product datasheet for **RR206495**

Hspa1l (NM_212546) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Hspa1l (NM_212546) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Hspa1l
Synonyms:	Hsp70-3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR206495 representing NM_212546
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCTGCTAATAAAGGAATGGCCATAGGCATCGACCTGGGCACCACCTACTCGTGGTGGCGTGTTC
 AGCACGGCAAGGTGGAGATCATCGCCAACGACCAGGGCAACCGCACGACCCCAAGCTACGTGGCCTTAC
 CGACACCGAGCGGCTCATCGGGGACGCCCAAGAACCAGGTGGCCATGAATCCCCAGAACACTGTTTTT
 GACGCCAAACGTCTAATTGGCAGGAAGTTAATGATCCCGTCGTGCAGTCGGATATGAAGCTGTGGCCAT
 TTCAAGTGATCAATGAAGCGGGCAAGCCTAAGGTGCTGGTGTCTACAAAGGAGAGAAGAAAGCCTTCTA
 CCCGGAGGAGATCTCGTCCATGGTGTGACCAAGATGAAGGAGACTGCAGAGGCTTTTTTGGCCACAGT
 GTCACCAACGCTGTGATCACCGTCCAGCCTATTTCAATGACTCCAGAGACAGGCCACTAAGATGCAG
 GGGTCATCGCAGGACTCAACGTGCTGCGAATCATCAATGAGCCCACAGCGGGCCATCGCTACGGCTT
 GGATAAAGGAAGTCACGGAGAGCGGCACGTGCTCATCTTCGACCTGGGGGCGGCACGTTTCGACGTGCC
 ATCCTGACGATCGACGACGGCATCTTCGAGGTGAAGGCCACGGCGGGCGACACGCACCTGGCGGGGAGG
 ACTTCGACAACCGGCTGGTGGCCACTTCGTGGAGGAGTTCAAGAGGAAGCACAAGAAGGACATCAGCCA
 GAACAAGCGCGCGGTGCGGCGTCTGCGCACGGCGTGCAGAGGGCCAAGAGGACGCTGTCGTCCAGCACC
 CAGGCCAACCTGGAGATCGACTCTCTGTACGAGGGCATCGACTTCTACACGTCCATCACGCGGGCGCGGT
 TCGAGGAGCTGTGTGCAGACCTATTTAGAGGCACGCTTGAGCCCGTGGAGAAGTCTCTCCGGGACGCCAA
 GATGGACAAGGCTAAAATCCATGACATTGTTCTAGTAGGGGGCTTACCCGCATCCCAAAGGTGCAAAAAG
 CTGCTTCAAGACTACTTTAATGGACGGGATCTCAACAAGAGTATCAATCCCGATGAGGCGGTGGCGTACG
 CGTAGCTCCCCTGTCTCTAGGATTGGAGACAGCCGGGGGTGTGATGACTGTTCTGATCAAGCGCAACTCC
 ACCATCCCACCAAGCAGACGCAGATCTTCACCACCTACTCGGACAACCAGCCCGGGGTGCTGATCCAGG
 TGTACGAGGGCGAGAGGGCCATGACGCGCGACAACAACCTCCTGGGGCGCTTTGACCTGACTGGAATACC
 CCCTGCACCTAGGGGTGTGCCCCAGATCGAGGTGACCTTCGACATCGACGCCAACGGGATTCTCAACGTC
 ACAGCCATGGACAAGAGCACCGGCAAGGCCAACAAGATCACCATCACCACGACAAGGGCCGCTTGAGCA
 AGGAGGAGATCGAGCGCATGGTGAAGAGGCCGAGCGGTACAAGGCCGAGGACGAGGGCCAGAGGGAGAA
 GATCGCTGCCAAAATGCCTTAGAATCGTATGCCTTTAACATGAAGAGCGCTGTAGGCGATGAGGGTCTG
 AAGGACAAGATCAGCGAGTCTGATAAAAAGAAAATACTGGATAAATGCAGTGAGGTCCTTCTCTGGCTGG
 AGGCCAACCGAGCTGGCGGAGAAAGAAGAGTTTCGATCATAAAAGAAAAGAGCTGGAAAACATGTGCAACCC
 TATCATACAAAGCTGTACCAGAGCGGCTGCACAGGGCCACCTGCGCGCCAGGGTATACACCCGGCAGG
 GCTGCCACAGGCCCCACCATTGAGGAAGTAGAT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206495 representing NM_212546
 Red=Cloning site Green=Tags(s)


MAANKGMAIGIDLGTTYSCVGVFQHGKVEIIANDQGNRTTPSYVAFDTERLIGDAAKNQVAMNPQNTVF
 DAKRLIGRKFNDPVVQSDMKLWPFQVINEAGPKVLVSYKGEKKAFYPPEISSMVLTKMKETAEAF LGHS
 VTNAVITVPAYFNDSQRQATKDAGVIAGLNLVRIINEPTAAAIAYGLDKGSHGERHVLIFDLGGGTFDVS
 ILTIDDGIFEVKATAGDTHLGGEDFDNRLVSHFVEEFKRKHKKDISQNKRAVRRRLTACERAKRTLSSST
 QANLEIDSLYEGIDFYTSITRARFEELCADLFRGTLEPVEKSLRDAKMDKAKIHDIIVLVGGSTRIPKVQK
 LLQDYFNDRDLNKSINPDEAVAYGAAVQAAIILMGDKSEKVQDLLLLDVAPLSLGLTAGGVMTVLIKRN
 TIPTKQTQIFTTYSDNQPGLIQVYEGERAMTRDNLLGRFDLTGIPPAPRGVPIEVTFFDIDANGILNV
 TAMDKSTGKANKITITNDKGRLSKEEIERMVQEAERYKADEGQREKIAAKNALESYAFNMKSAVGEGL
 KDKISESDKKKILDKCSEVLSWLEANQLAEKEEFDHKRKELENMCNPIITKLYQSGCTGPTCAPGYTPGR
 AATGPTIEEVD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_212546

ORF Size: 1923 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_212546.4](#), [NP_997711.1](#)

RefSeq Size:	2632 bp
RefSeq ORF:	1926 bp
Locus ID:	24963
UniProt ID:	P55063
Cytogenetics:	20p12
MW:	70.5 kDa
Gene Summary:	one of three major histocompatibility complex (MHC)-linked heat shock protein 70 (Hsp70) genes [RGD, Feb 2006]