

Product datasheet for **RR206527**

Rdh13 (NM_001108468) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rdh13 (NM_001108468) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Rdh13
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR206527 representing NM_001108468 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAATCGCTTTCTGCTGCCGGTGTCCGTAGTGGGCACGGTAATCGGCGGCACCGTGTCTCAAGGACT
ATGTCGCTGGTGGGCTTGTCCCAGCAAGGCCACCATAACCTGGGAGGACTGTTATCGTGACAGGAGCCAA
CACAGGCATCGGAAACAACTGCTTTGGAGCTGGCTAAAAGAGGAGGCAACATCATTCTGCGCTGCCGT
GACAGGGAGAAGTGTGAGGCAGCTGCAAAGGACATTCGAGGAGAGACCTTGAATCCCGTGTGCGTGCCG
AGCACCTAGACTTGGCTTCCCTCAAGTCCATTCGAGAGTTTGCAGGGAAGATCATTAAAGAGGAGGAGCG
AGTAGACATTCTGGTCAACAATGCAGCAGTGATGCGGTGCCCACTGGACCACCGAGGATGGCTTTGAG
ATGCAGTTTGGTGTCAACCACTTGGGCAATTTCTATTGACAAATTTGCTGTTGGACAAGCTGAAGGCCCT
CGGCCCATCACGCATCATCAATCTCCTCCCTGGCTCATGTTGCTGGGCACATAGACTTCGAGGACTT
GAACTGGCAGATAAAAAAGTATGACACCAAGGCAGCTTACTGCCAGAGCAAGTTGGCTGTTGCTCCTCTTC
ACCAAGGAGCTGAGTCGCCGGCTGCAAGGCACTGGTGTGACTGTCAATGCCTGCACCCTGGCGTGGCCA
GGACGGAGCTGGGGCAGACATACAGGCATGCACAACTCTGCGTTCTCTGGCTTCATGCTTGGGCCCTCTT
CTGGCTGCTGTTCAAGAGTCCCCAGCTGGCGGCCAGCCAGCACATAACCTGGCTGTGGCAGAGGAACTG
GAGAGTGTCTCTGGGAAGTACTTTGATGGACTCAGAGAGAAGGCTCCATCTCCTGAGGCTGAAGATGAGG
AAGTAGCCCGGAGGCTTTGGACTGAAAGTGCCATTTGGTGGGCTTGGACATGGCTCATGGGTCCTCTGG
GAGAGGACATTCATCTCCAGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR206527 representing NM_001108468
 Red=Cloning site Green=Tags(s)

MNRFLLPVSVVGTIVIGTVLLKDYVAGGACPSKATIPGRTVIVTGANTGIGKQTALELAKRGGNIILACR
 DREKCEAAAKDIRGETLNPRVRAEHLDLASLKSIREFAGKIIKEEERVDILVNNAVMRCPHWTTEDGFE
 MQFGVNHLDGHFLLTNLLLDKLGKASAPSRRIINLSSLAHVAGHIDFDLNLWQIKKYDTKAAQCQSKLAVVLF
 TKELSRRLQGTGVTVNALHPGVARTELGRHTGMHNSAFSGFMLGPFFWLLFKSPQLAAQPSTYLVAEEL
 ESVSGKYFDGLREKAPSPEAEDEEVARRLWTESAHLVGLDMAHGSSGRGHSISR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108468

ORF Size:	1002 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001108468.1 , NP_001101938.1
RefSeq Size:	2632 bp
RefSeq ORF:	1005 bp
Locus ID:	361504
Cytogenetics:	1q12
MW:	36.4 kDa