

Product datasheet for **RR206542L3V**

Nedd4l (NM_001008300) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Nedd4l (NM_001008300) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Nedd4l
Synonyms:	Nedd4-2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001008300
ORF Size:	2889 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206542).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008300.2 , NP_001008301.2
RefSeq Size:	4213 bp
RefSeq ORF:	2892 bp
Locus ID:	291553
Cytogenetics:	18q12.1
Gene Summary:	mouse homolog appears to be involved in binding and regulating the amiloride-sensitive epithelial sodium channel which plays a role in sodium homeostasis and control of blood volume and blood pressure [RGD, Feb 2006]



[View online »](#)