

Product datasheet for **RR206598**

Slc15a3 (NM_139341) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc15a3 (NM_139341) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc15a3
Synonyms:	Pht2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>RR206598 representing NM_139341
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGTCGCACTTCGCGCCGAACAGCAGCAAGCAGATCCGGAGAGAGGCAACCACTGGTGGCTCAAGGCC
 GGTGGGGGCCCGACGGTGGAGGCGGACAGCAGCTGCTGCGGTGCTGCTAGTGGAGATGCTGGAACGCGC
 TGCCTTCTTTGGCGTCACTCCAACCTCGTGCTCTACCTCAACAGCTTGAACCTCAACTGGGATGGTGAG
 CACGCGTCGCGGGCCACACTCCTCTTCTTAGGGCCTCCTACCTGCTGGCTCCGGTGGGAGGCTGGCTGG
 CCGACGTGTACCTGGGCCGTTTCTTGTATTTTCGCTCAGCCTGCTGCTCTACCTGGCTGCCACTGGCTT
 GTTGCTCACCACCATCACCGACGATGGCCGAGATCCTTCTGTGGAGAGATGCCCGAGTTGCCACTGAAA
 CCTGCCTGTCCATCCGCAAACCTGCAAGGCTCCTGGTCCAGCCCCTACTGCGCGACCACCCTCTACTTAG
 TGTGCTTCTCTTGGCGCTGGTCCAGCTCCGTGAGGAGCAACCTCACCTCGTTCGGGGCCGACCAGGT
 GATGGATCTCGGCCGTGATGCCACCCGCGCTTCTTCAACTGGTCTATTGGAGCATCAATCTGGGTGCC
 ATTCTGTCCCTGCTGGTGGTGGCTTCAATCGAACAGAATCAGCTTCTACAGGGCTACAGCATCATCG
 TGGGCCTCGTGGGCTGGCATTCTCATTCTTCTCATCGCCACACCCGCTTTCATCACCAGCCCCAAC
 AGGCAGCAAGTGTATCTATGCTGAACCTTGCAATCCAAAACCTGCTGTCCCGGGTGGCAGTGGTGGCGG
 AGGCCCTCTCCAGGAACCTGAAGGTGCTCACCTGTTGCCTGACCAGAGGTCTAACCAGCCTGGACCTT
 CTCCACAAGAAGATAGGCCAACTTCCAGGTGCTGTGAAGGTCTGCCTGTGATGGTGACCCCTGTGCTC
 TACTGGATGGTGTATTTCCAGATGCAGTCCACCTATGCTCTACAAGGTCTCCACCTCCATATCCCCAAC
 ATCTCCGCACCAACCCTAACATCTCTTTGCCACTGAGATCGGATAGCAGTAACCTACAGGATCCAGAAG
 CCTGGCTCCTGCTGGCAATGTTGCAGTGATACTGATTCTGGTCCCTGTCAAGGACCCTAATTTGATCC
 TCTGCTGCTGCGGTGTAAGCTGCTGCCCTCGGCTCTGCAGAAGATGGCGCTGGGATGTTCTTTGGGTTT
 ACTTCCATCATCGTGGCAGGAGTCTGAAAAAGAGCGCTTGCAGTACATTGCTGCCAACAGACAGTGC
 CTCAGCTGATTGGGAAGGACTTATACTATGCGGCACCACTGTCCATCTGGTGGCAAATCCCCAATACCT
 GCTCATCGGGATCAGTGAGATTTTCCAGCATCCAGGTCTGGAGTTCGCTTACTCGGAAGCCCCACGC
 TCCATGCAGGGCGCAATCATGGGCATCTTCTTGCCTTTCGGGGGTAGGTTCACTGCTAGGCTCTGGCT
 TGGTGGCCCTACTGTCCCTTCTGGAGGCTGGATGTACTGCCCTAAGGACTTCGGGAACATCAACAACCTG
 TCGTATGGACCTTACTTCTTCTGCTAGCCGGCATTGAGGCTGTACAGCTGTCTGTTCTCTGGATT
 GCTGGTGTATGAGAGGACTCGTCAGGACCCAGACTCCAGAACTCAACCAGCAGAGTCAAGGGC

AG**CGGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>RR206598 representing NM_139341
 Red=Cloning site Green=Tags(s)


MSALRAEQPSRSGERQPLVAQGRWGP RRWRRTAAAAVLLVEMLERAFFGVTSNLVLYLNSLNFNWDGE
 HASRATLLFLGASYLLAPVGGWLADVYLGRFLAISL SLLL YLAATGLLLTTITDDGRRSFCGEMPELPLK
 PACPSANCQGSWSSPYCATTLYL VLLLLALAASSVRSNLT SFGADQVMDLGRDATRRFFNWFYWSINLGA
 ILSLLVVAFIEQNISFLQGSIIIVGLVGLAFFIFLIATPVFITKPP TGSQVSSMLNLA FNCCPGWQWWR
 RPSSRNSEGAHLLPDQRSNQPGSPQEDMANFQVLLKVL PVMVTLVPYWMVYFQMSTYVLQGLHLHIPN
 IFRTNPNISLPLRSDSSNYRIPEAWLLANVAVILILVPVKDHLIDPLLLRCKLLPSALQKMG MFFGF
 TSIIVAGVLEKERLQYIAANQTPQLIGKDLYYAAPLSIWWQIPQYLLIGISEIFASIPGLEFAYSEAPR
 SMQGAIMGIFFLSGVGSLLGSGLVALLSLPGGWMYCPKDFGNINCRMDLYFFLLAGIQAVTAVLFLWI
 AGRYERTRQDPDSQNSTSRVRG

SGP**TRTRPLEQKLI**SEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_139341
 ORF Size: 1746 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_139341.1 , NP_647557.1
RefSeq Size:	1979 bp
RefSeq ORF:	1749 bp
Locus ID:	246239
UniProt ID:	Q924V4
Cytogenetics:	1q43
MW:	64.6 kDa
Gene Summary:	may play a role in protein catabolism in the lymphatic system [RGD, Feb 2006]