

## Product datasheet for **RR206605L4V**

### **Kcne3 (NM\_022235) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Kcne3 (NM_022235) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Kcne3
Synonyms:	MGC105432
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022235
ORF Size:	321 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206605).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022235.2</a> , <a href="#">NP_071571.1</a>
RefSeq Size:	1205 bp
RefSeq ORF:	324 bp
Locus ID:	63883
UniProt ID:	<a href="#">Q9JIV7</a>
Cytogenetics:	1q32
Gene Summary:	a regulator that is coexpressed with K+ channel KCNQ1 [RGD, Feb 2006]



[View online »](#)