

Product datasheet for **RR206607L4V**

Eif1ad (NM_001008305) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Eif1ad (NM_001008305) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Eif1ad
Synonyms:	RGD1304686
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008305
ORF Size:	501 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206607).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008305.1 , NP_001008306.1
RefSeq Size:	2577 bp
RefSeq ORF:	504 bp
Locus ID:	293673
UniProt ID:	Q5RKI6
Cytogenetics:	1q43
Gene Summary:	Plays a role into cellular response to oxidative stress. Decreases cell proliferation (By similarity).[UniProtKB/Swiss-Prot Function]



[View online »](#)