

Product datasheet for RR206623

Noxo1 (NM_001106986) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Noxo1 (NM_001106986) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Noxo1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206623 representing NM_001106986
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGCAAGTCCAAGACACCCAGTATCAGCCCATGCTGTAGCCCTGGTGCAGATGGAACGGCTCCAGACGT
TTGCCTTCTCTGTGTGGTCAGACAACAGTGACACCTTTGCTCGCAGGAGCTGGGAGGAGTTTAGCA
GCTCCAGAAGACCCTTAAGAAAATCTCCAGTGGAGGCAGGCCCTGCTTCAGAGATCTGAACGAGTCTT
CCCAAGCTTCTGATGCTCCATTGTTGACACGTCGGGGACATACTGGCCGAGCCCTGTTACGTTTGC
TGTGGAAACCTATGTACGGTCACTGCTGGCAACCTCACAGCACATAGTGACGAGCTCAACACTCAATAG
CTTCTTTGCACCAAACTCTGGATCTGGAGCCCATGCTGCCTCCTGGCAGCCTGGTATCCTGCCACC
CCAGAGGAGCCCTTATCCAACCTATAGGCAGCCTTGCCATTATAGCCTGGAGGCTCAGAGCATGCGCT
GTTACAGCCATTCCACTCTTGACACAAAAGACAGACCTTTCCACTCAGGCTAAAGAAATTCTGGA
CATATTACTACGGCATCCTTCAGGCTGGTGGTGGTGGAGAACAAGGACCAGCAGACAGCCTGGTTTCCA
GCTCCCTACCTGGAGGAGATAGCAACAGGCCAAGGCCAGGAGTCGGGCATGGCTGTTCAAGGAAGTGTA
TGAGC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206623 representing NM_001106986
 Red=Cloning site Green=Tags(s)

MASPRHPVSAHAVALVQMERLQTFAFSVCWSDNSDTFARRSWEEFRQLQKTLKKIFPVEAGLLQRSERVL
PKLPDAPLLTRRGHTGRGLLRLLRLETYVRSLLATSQHIYTSSTLNSFFAPKPLDLEPMLPPGSLVILPT
PEEPLSQPIGSLAIHSLEAQMRLQPFHTLDTKDRPFHTQAKEILDILLRHPSGWWLVENKDDQQTAWFP
APYLEEIAATGQQQESGMAVQGGMS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV




[View online »](#)

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106986

ORF Size: 705 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001106986.1, NP_001100456.1</u>
RefSeq Size:	1965 bp
RefSeq ORF:	708 bp
Locus ID:	302976
Cytogenetics:	10q12
MW:	26.4 kDa