

Product datasheet for **RR206630**

Ogn (NM_001106103) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Ogn (NM_001106103) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Ogn
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR206630 representing NM_001106103
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGACTGTGCATCCTACGCTTCTCCTGCTATTGTTTCGTGCCTCTGACACAGCAAACACCACAAACTC
AGCTGGACTCACACGTTAACTCTAAGTATGGAACAGATAATTCTGAAGAAACAAATTTATCCAGGATTA
TGAAGATAAATACCTGGATGGGAAAAATATTAAGGAAAAAGAACTATGATAATTCCTGATGAGAAAAGT
CTTCAGTTACAAAAGATGAAGTCATACCATCAGTGCCACCAAGAAAGAAAATGATGAAATGCCACAT
GCCTGTTGTGTCTGTTAAGTGGCTCTGTCTACTGTGAAGAAGTTGACATTGATGCTGTACCACCACT
GCCAAAAGAATCAGCTTATCTTTATGCACGATTCAACAAAATAAAAAGCTGACTGCCAAGGATTTTGCA
GACATGCCTAACTTAAGAAGACTTGATTTTACGGGAAATTTGATAGAAGACATAGAAGATGGTACTTTTT
CAAACTTTCTCTGTTAGAAGAACTTACACTTGCTGAAAACCACTACTAAGACTTCCAGTTCTTCTCC
AAAGCTTACTTTACTCAATGCCAAATACAACAAAATCAAGAGTAAAGGAATTAAGCAAACACATTCAA
AAACTGAATAAACTCTCTTTCTCTATTTGGACCATAATGATCTGGAATCTGTGCCTCCTAATTTACCAG
AAAGTCTACGTGAATTCACCTTCAGTTAACAGCATATCTTCAATTACCGATGATACATTCTGCAAGGC
TAATGACTCGGTACATTCGGGAGCGCATGGAAGAGATCCGCCTGGAGGCAATCCGATTGCTCTGGGA
AAGCATCAAACAGCTTTATCTGCTTAAAAAGATTACCCACAGGGTCATACCTT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR206630 representing NM_001106103
 Red=Cloning site Green=Tags(s)

MKTVHPTLLLLLVPLTQQTPQTQLDSHVNSKYGTDNSEETKFIQDYEDKYL DGKNIKEKETMIIPDEKS
 LQLQKDEVIPSVPTKKENDEMP TCLLCVCLSGSVYCEEVDIDAVPPLPKESAYLYARFNKIKKLTAKDFA
 DMPNLRRLDFTGNLI EDIEDGTFSKLSLLEELTLAENQLRLPVLPPKLLNNAKYNKIKSKGIKANTFK
 KLNKLSFLYLDHNDLESVPPNLPESLRVIHLQFNSSITDDTFCKANDTRYIRERMEEIRLEGNPIALG
 KHPNSFICLKRLLPTGSYF


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001106103

ORF Size:	894 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001106103.1 , NP_001099573.1
RefSeq Size:	2650 bp
RefSeq ORF:	897 bp
Locus ID:	291015
Cytogenetics:	17p14
MW:	34.1 kDa