

Product datasheet for **RR206646L3V**

Cdc37I1 (NM_001011941) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdc37I1 (NM_001011941) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cdc37I1
Synonyms:	Cdc37I
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001011941
ORF Size:	1005 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206646).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001011941.1 , NP_001011941.1
RefSeq Size:	1633 bp
RefSeq ORF:	1008 bp
Locus ID:	293886
UniProt ID:	Q5XIC3
Cytogenetics:	1q52
Gene Summary:	Co-chaperone that binds to numerous proteins and promotes their interaction with Hsp70 and Hsp90.[UniProtKB/Swiss-Prot Function]



[View online »](#)