

Product datasheet for **RR206650L4V**

Grb2 (NM_030846) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Grb2 (NM_030846) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Grb2
Synonyms:	Ash-psi
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030846
ORF Size:	651 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206650).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030846.2 , NP_110473.2
RefSeq Size:	2608 bp
RefSeq ORF:	654 bp
Locus ID:	81504
UniProt ID:	P62994
Cytogenetics:	10q32.1
Gene Summary:	acts as a link between tyrosine kinase receptors and Ras signaling [RGD, Feb 2006]



[View online »](#)