

Product datasheet for **RR206662L4V**

Cdkn2c (NM_131902) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cdkn2c (NM_131902) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cdkn2c
Synonyms:	p18
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_131902
ORF Size:	504 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206662).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_131902.3 , NP_571977.1
RefSeq Size:	1462 bp
RefSeq ORF:	507 bp
Locus ID:	54238
Cytogenetics:	5q35
Gene Summary:	cyclin-dependent kinase inhibitor that may be involved in cell cycle regulation during terminal differentiation of oligodendrocytes [RGD, Feb 2006]



[View online »](#)