

Product datasheet for **RR206670L3V**

Dcps (NM_153302) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dcps (NM_153302) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dcps
Synonyms:	Hint-5
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_153302
ORF Size:	1008 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206670).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_153302.1 , NP_695214.1
RefSeq Size:	1221 bp
RefSeq ORF:	1011 bp
Locus ID:	266605
UniProt ID:	Q8K4F7
Cytogenetics:	8q21
Gene Summary:	human homolog is a scavenger pyrophosphatase that hydrolyzes the residual cap structure following 3 to 5 decay of an mRNA [RGD, Feb 2006]



[View online »](#)