

Product datasheet for **RR206711L4V**

Dab2ip (NM_138710) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Dab2ip (NM_138710) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dab2ip
Synonyms:	AIP-1; DIP1/2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138710
ORF Size:	3483 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206711).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138710.3 , NP_619724.3
RefSeq Size:	6508 bp
RefSeq ORF:	3486 bp
Locus ID:	192126
UniProt ID:	Q6P730
Cytogenetics:	3p11
Gene Summary:	GTPase-activating protein; binds the N-terminal domain of DOC-2/DAB2 [RGD, Feb 2006]



[View online »](#)