

## Product datasheet for **RR206714**

### Rab4b (NM\_017355) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Rab4b (NM\_017355) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Rab4b  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR206714 representing NM\_017355  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGGCTGAGACCTACGACTTCTCTTCAAATTCCTGGTGATTGGCAGTGCCGGAAGTGGCAATCATGTC  
TCCTCCATCAGTTTATTGAGAATAAGTTCAAACAGGACTCAACCCACACTATTGGCGTGGAGTTTGGATC  
CAGGGTTGTCAACGTTGGGGGAAGACTGTGAAGCTACAGATTTGGGATACTGCTGGCCAGGAGCGGTTT  
CGGTCCGTGACACGGAGTTATTACCGAGGGCAGCTGGAGCCCTGCTGGTGTATGACATACCAGCCGGG  
AGACATACAACTCGCTAGCTGCTTGGCTGACCGATGCCCGGACGCTGGCCAGCCCAACATCGTGGTCAT  
CCTCTGTGGCAACAAGAAGGACCTGGATCCTGAGCGTGAGGTCACTTTCTAGAGGCTTCCCGCTTGGCT  
CAGGAGAATGAGCTGATGTTCTGGAGACCAGCGCTCTCACGGGTGAAAACGTGGAAGAAGCTTTCTGA  
AGTGCGCACGCACCATCCTCAACAAGATCGACTCAGGTGAACTAGATCCGGAGAGGATGGGCTCAGGCAT  
TCAGTATGGTGACATATCTTTCGCCAGCTGCGGCAGCCTCGCAGTGCCAGGCCGTGGCCCTCAGCC  
TGTGGCTGC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >RR206714 representing NM\_017355  
Red=Cloning site Green=Tags(s)

MAETYDFLFKFLVIGSAGTGKSCLLHQFIENKFKQDSNHTIGVEFGSRVNVVGGKTVKLIWDTAGQERF  
RSVTRSYYRGAAGALLYVDITSRETYNSLAAWLTDARTLSPNIVVILCGNKKDLDPEREVTFLEASRFA  
QENELMFLSALTGENVEEAFLKCARTILNKIDSSELDPERMGSGIQYGDISLRQLRQPRSAQAVAPQP  
CGC

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV



[View online »](#)



<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_017355.1</a></u> , <u><a href="#">NP_059051.1</a></u>
<b>RefSeq Size:</b>	1130 bp
<b>RefSeq ORF:</b>	642 bp
<b>Locus ID:</b>	50866
<b>UniProt ID:</b>	<u><a href="#">P51146</a></u>
<b>Cytogenetics:</b>	1q21
<b>MW:</b>	23.6 kDa
<b>Gene Summary:</b>	may play a role in regulation of membrane trafficking [RGD, Feb 2006]