

Product datasheet for **RR206736L4V**

Alas2 (NM_013197) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Alas2 (NM_013197) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Alas2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_013197
ORF Size:	1761 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206736).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_013197.1 , NP_037329.1
RefSeq Size:	1899 bp
RefSeq ORF:	1764 bp
Locus ID:	25748
UniProt ID:	Q63147
Cytogenetics:	Xq13
Gene Summary:	erythroid-specific isozyme that catalyzes the conversion of glycine and succinyl-CoA to 5-aminolevulinate [RGD, Feb 2006]



[View online »](#)