

Product datasheet for **RR206749L3V**

Fut10 (NM_001012360) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fut10 (NM_001012360) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fut10
Synonyms:	fut10A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001012360
ORF Size:	1308 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206749).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001012360.1 , NP_001012360.1
RefSeq Size:	1311 bp
RefSeq ORF:	1311 bp
Locus ID:	497619
UniProt ID:	Q5F2L1
Cytogenetics:	16q12.3
Gene Summary:	Probable fucosyltransferase.[UniProtKB/Swiss-Prot Function]



[View online »](#)