

Product datasheet for **RR206752L4V**

Cd93 (NM_053383) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cd93 (NM_053383) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cd93
Synonyms:	C1qr1; C1qRp; Ly68
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053383
ORF Size:	1929 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206752).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053383.1
RefSeq Size:	4779 bp
RefSeq ORF:	1932 bp
Locus ID:	84398
UniProt ID:	Q9ET61
Cytogenetics:	3q41
Gene Summary:	enhances cell-cell interactions; may play a role in innate immune response against pathogens and toxic cell debris [RGD, Feb 2006]



[View online »](#)