

Product datasheet for **RR206759**

Slc6a14 (NM_001037544) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Slc6a14 (NM_001037544) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Slc6a14
Synonyms:	Atb0+
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>RR206759 representing NM_001037544
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGACAGATTGAAGTGCCCGAGCTTCTCAAGTGCAGACAGAAGGAGAAAAGTGACAGCTTCATCTGAGA
ACTTCCATGTTGGTGAAAATGATGAAAATCAGGAGCGTGGTAACTGGTCCAAAAAATCAGATTATCTTCT
TCCATGTTGGCTATGCAGTGGGATTGGGCAACGTGTGGCGATTCCCTTACCTGACCTACACCAATGGT
GGAGGTGCTTTCTTGATACCCTATGCAATTATGCTAGCATTGGCTGGCTTGCCTTTGTTCTTCTTGAAT
GTTCACTGGGACAATTTGCTAGCTTAGGTCCAGTTTCAGTTTGGAGGATTTTGCCTGTTTCAAGGTGT
GGGAATCACGATGGTCTGATATCCATTTTTGTGACAATTTATTATAATGCATCATTGCCTATAGTCTT
TACTACTATTTGCTTCTTTCAAAGTGTACTGCCATGGGCTAATTGTTCTTCTTGGCAGACGATA
ACTGTAGCAGACAACCAATAGTAACTGGCTGTCGTGTGAGCACAGACACAGGAGGAGACATGTTTATGAA
TATAAGCTGGGTAATAGCAACAACCTCACCTGCTTGAATGGTAGTGAAGTTTTTCGGCCAGGACAACCT
CCCAGTGAACAGTATTGGGATAAAGTGACGCTTCGGAGGTCAAGTGGAATGGATGAAACTGGCGTAATTG
TGTGGTACTTAGCACTTTGCCTTCTTAGCTGGATCATAGTTGGAGCAGCACTATTTAAAGGAATCAA
ATCCTCTGGCAAGGTGGTATATTTACAGCTCTTTCCCGTATGTGGTCTTCTCATTCTACTAATACGA
GGTGCACTTTGGAGGTGCATCAAAGGCATTTCATATTACATTGGAGCACAGTCAAATTTACAAAAC
TCAGAGAAGCTGAGGTTTGAAGGATGCTGCCACTCAGATATTTTATCCCTTTCTGTGGCTTGGGGTGG
CCTAGTGGCTCTGTCATTTACAACAAGTTCAACAACAACCTGTTACTCTGATGCCATTATAGTTTGTCTG
ACCAACTGCCTCACTAGTGTGTTTGTGGATTGCGCATTTTTTCTATTCTGGACACATGGCTCATATAT
CCGGAAGGAAGTCTCTCAAGTTGAAAAATCAGGATTTGACTTGGCATTTCATTGCTTATCCGGAAGCACT
AGCTCAACTCCCAGGTGGTCCATTCTGGTCCATATTATTTTTCTCATGCTTTTAAACATTGGGTCTGGAT
TCTCAATTTGCTTCTATTGAAACAATCACAACAACAATCAAGATTTGTTTCCAAAGCAATGAAACGCA
TGAGGGTCTCTATAACTTTGGGATGCTGCTTGTATTTTGTCTCCTTGGTCTCGTCTGTGTGACTCAGGC
TGGAAATTTACTGGGTTCAATTTGATTGATCACTTCTGTGCTGGATGGGGATCTTAATTTGCTGCCATCTT
GAAATTCAGGAATCATCTGGATTTATGGAGGAAATAGGTTTATTGAAGATATAGAAATGATGATTGGAG
CAAAAAGATGGCTATTCTGGCTTTGGTGGAGAGCTTGTGGTTTGTGTCATCACTCCTATACTGTTGTCTGC
AATTCTAGTCTGGTCACTGGTGAATTTACAGACCTAATTATGCTGATATCCATATCCTGACTGGGGA
GTCGCTCTCGGATGGTGTATGATTATTTCTGCATTATTTGGATACCAATTATGGCGATCATAAAAATAG
TTCGGGCTGAAGGAAACATCCTTCAGCGTATCATAAAATGCTGTAGACCAGCTTCAAACCTGGGGTCCATA
CCTGGA AAAACATCGTGGAGAGAGATATAAAGACATGGCTGAACCTGCAAAGGAGACTGACCACGAAATA
CCTACTATTAGTGGCACTAGAAAACCAGAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR206759 representing NM_001037544
Red=Cloning site Green=Tags(s)


```
MDRLKCPSPFFKCRQKEKVTASSENFHVGENDENQERGNWSKSDYLLSMVGYAVGLGNVWRFPPLYTYTNG
GGAFILIPYAIMLALAGLPLFFLECSLGFASLGPVSVWRILPLFQGVGITMVLISIFVTIYYNVIAYS
YYLFAFQSVLPWANCSSWADDNCSRQPIVTGCRVSTDTGGDMFMNISWVNSNNFTCLNGSEVFRPQQL
PSEQYWDKVTLLRRSSGMDGTGVIWYLALCLLLAWIIVGAALFKGIKSSGKVYVFTALFPYVLLILLIR
GATLEGASKGISYYIGAQSNFTKLREAEVWKDAATQIFYSLVAVWGLVALSSYNKFNNNCYSDAIIIVCL
TNCLTSVFAGFAIFSILGHMAHISGKEVSQVVKSGFDLAFIAYPEALQPLGGPFWSILFFFMLLLGLD
SQFASIIETITTTIQDLFPKAMKMRVPITLGCCLILFLLGLVCVTQAGIYWVHLIDHFCAGWGILIAAIL
EIAGIIWIYGGNRFIEDIEMMIGAKRWLFWLWWRACWFVITPILLSAILVWLVKFRPNYADIPYPDWG
VALGWCMIIFCIWIPIMAIKIVRAEGNILQRIKCCRPASNWGPYLEKHRGERYKDMAEPAKETDHEI
PTISGTRKPE
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001037544

ORF Size: 1920 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001037544.1](#), [NP_001032633.1](#)

RefSeq Size: 2138 bp


RefSeq ORF: 1923 bp

Locus ID: 298340

Cytogenetics: Xq34

MW: 71.8 kDa

Product images:



Circular map for RR206759