

## Product datasheet for **RR206834**

### Fam234a (NM\_001009701) Rat Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Fam234a (NM_001009701) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Fam234a
Synonyms:	Itfg3; RGD1311484
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

**ORF Nucleotide Sequence:**

>RR206834 representing NM\_001009701  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGATGGATGACAAGGATCTAGAAGCTGAAATACATCCCTTGAAGAATGAGGACAAAAAGTCACAGGAAA  
 ACCTAGGAAACCTACCAAGACCGAGGACAACCTGAAGAACAAGCCTGTGCCTTCTCGCTTGTCCCGGTG  
 TCGCACGGTGGCATTTCCTTTCCCTGTTTCATCTGTCTCTTTGTGGTGTGGTCTCTTTTCATCATC  
 CCATGCCAGACCGACCGTCTCTGAGGACAAGTGGAGGCTGGACTACGACTATGCAGTCACGTACGACT  
 TTCTGGCTCTGGGAGACGTAGACAGGGACGAGGTTGAGGACGTTCTTTCTTTATAAAAACACCAACAG  
 CAGTAACAATCCACCAGGTCTGTGCTGACGAAGGCTTTTCCACTCCCTGCACCTTTATGGCTGCTGTG  
 TCAGGGACCAACGGTACCGTCTCTGGGAGCGCCTTTGGCCCAAGATGTTGCCCATGTAAGTGTGCCG  
 TGCCACAGACACGGGACAGCAAGGTGTCTTCTGCCTGTGTCTGGTGGGAAGACTGGTTCTTTTCATGGC  
 AGTCAACTTCTTACAGGAGAAACGCTGTGGAGCCATCCAGAAGCTTTAGTGGGAATGCGTCCATCTTG  
 AGTCTCTGCTCCAGGTGCCTGACATCGACGGTACGGTGTCCAGACCTGCTGATGCTCTCCCGGGAGG  
 GACAGGAGGTCACTGGCGCTCTATTACAGGTAGCACTGGTACCAGATTGGCCACAGAGGCAGCCTTGG  
 TGTGGATGGAGATGGTGTGCGCCCTCTTACAGTACCAGAACGGGTGCTCAGTATATTCTCTGCGCTGT  
 GCAAGTGCCTCTGTGGCGTCTCTGTGAAGGTCTCTATAAGAGAATTACTGGGAGAGATGACCATTTTA  
 AGGAAGACCCCTCTGGGAGAACATGCTCAATCGCTCTGCTCACAGGAGGCTTTTGCACAGCTCTGGAGC  
 AGTTCGATACCTGATGAATGTCCAGGGAAAACGGCCAGGACTTGCTTCTTGTGAGCTCTGAGGACTGT  
 AAGCTGCTGGATGGACAAGATCTGGTCCCAGGTGGACTCTTGGTGAACCTCAGGTTCTGAGAAAACCTA  
 TCCTTGGCCACTACAAACCAGACACCTTGGCTGTGGTAATTGAAAATGGAACAGCCTTGTAGACAACTA  
 CCTCCTTCTGGACCTCAGCACTGGGTCTATGCTGTGGAGCCAGCCCTCCAAGCCTTCTGGGGCCCCA  
 CCATCCACAAGCCTGATGACTGCAGACCACCGCTCAGCCTTCTTCTTCTGGGGCCTCCATGACCCTGTGA  
 GCGCTAATGAGATGGATCCCCAGGACACGCAGCACAGCCTGTATATGTTCCACCCGACCCTGCCCGCAT  
 CCTGCTGGAGCTAGCCAACGTGTCCGCCAACATCGTGGCCTTTGATGCGGTCTGTTTGAGCCAAGCCGC  
 CACGCTGCCTATGTGCTCCTGACGGGCCAGCAAGTTCAGATGTGCTGGCCTGGTCTCTCTGACCAAAC  
 ACAGAGTGCAGACCTTGTCCCAGGCAGTCGAGTATCCACTTGAGCGAAGGCAGCTCCGACAGTGACCA  
 GGCTATCCGGGACCGGTTCTCCCGACTACGGTACCGGAGTGAGGTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:**

>RR206834 representing NM\_001009701  
 Red=Cloning site Green=Tags(s)


MMDDKLEAEIHPLKNEDKKSQENLGNLPKTEDNLKNKPVPSRLSRCRTVAFFLSLFI CLFVVFVLSFII  
 PCPDRPSSSEDKWRLDYDYAVTYDFLALGDVDRDEVQDVLFLYKNTNSSNNSTRSCADEGFSTPCTFMAAV  
 SGTNGTVLWERPLAQDVAHVKCAVPQTRDSKVSSACVLVGRLGSMFMAVNF FTGETLWSHPRFSGNASIL  
 SPLLQVPDIDGDGAPDLLML SREGQEVSGALYSGSTGYQIGHRGS LGVDGDGVALLHVTRTGAQYILLPC  
 ASALCGVSVKGLYKRITGRDDHFKEDPSWENMLNRSARHRLHSSGAVRYLMNVPKGTGQDLLLVSEDC  
 KLLDGDQLVPRWTLGATQVLRKPI LGHYKPD TLAVVIENG TSLDRQILLDLSTGSMLWSQPLPSLPGAP  
 PSTSLMTADHRSAFFFWGLHDPV SANEMDPQDTQHSLYMFHPTLP GILLELANV SANIVAFDAVLFEPSR  
 HAAYVLLTGPASSDVPGLVSLTKHRVRDLVPGSRVIHLSEGSSDSQAI RDRFSRLRYRSEV

**TR**TRPLEQKLI SEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN:

NM\_001009701

<b>ORF Size:</b>	1656 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001009701.1</a> , <a href="#">NP_001009701.1</a>
<b>RefSeq Size:</b>	2652 bp
<b>RefSeq ORF:</b>	1659 bp
<b>Locus ID:</b>	360502
<b>UniProt ID:</b>	<a href="#">Q5M7W6</a>
<b>Cytogenetics:</b>	10q12
<b>MW:</b>	60.6 kDa