

Product datasheet for **RR206835L4V**

Fgl1 (NM_172010) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Fgl1 (NM_172010) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Fgl1
Synonyms:	Frep1; Lfire1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_172010
ORF Size:	801 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206835).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_172010.1 , NP_742007.1
RefSeq Size:	1141 bp
RefSeq ORF:	804 bp
Locus ID:	246186
Cytogenetics:	16q12.1
Gene Summary:	member of the fibrinogen family; appears to lack the platelet-binding and thrombin-sensitive sites, and cross-linking region, which are necessary for fibrin clot formation; protein may be liver-specific [RGD, Feb 2006]



[View online »](#)