

Product datasheet for **RR206837**

Stx1a (NM_053788) Rat Tagged ORF Clone

Product data:

| | |
|--------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Stx1a (NM_053788) Rat Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Stx1a |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| ORF Nucleotide Sequence: | >RR206837 representing NM_053788 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGAAGGACCGAACCAGGAGCTCCGCACGGCCAAGGACAGCGATGACGACGATGATGTCCTGTCCTG
TGGACCGAGACCGCTTCATGGATGAGTTCTTTGAACAGGTGGAAGAGATCCGGGGCTTTATTGACAAGAT
TGCTGAAAACGTGGAGGAGGTGAAGAGGAAACACAGCGCCATCCTGGCCTCCCCGAACCCGATGAGAAG
ACCAAGGAGGAACTGGAGGAGCTCATGTCGGACATTAAGAAGACAGCGAACAAGTTCGCTCCAAGCTAA
AGAGCATCGAGCAGAGCATCGAGCAGGAGGAAGTCTGAACCGCTCGTCGGCGGACCTGAGGATCCGGAA
GACGCAGCATTCCACGCTGTCCCGAAAAGTTGTGGAGGTCATGTCCGAGTACAACGCCACTCAGTCAGAC
TACCGAGAACGCTGCAAAGGGCGCATCCAGAGGCAGCTGGAGATCACTGGCCGGACCAGCAGCAGTGGG
AGTTGGAAGACATGCTGGAGAGTGGGAATCCCGCCATCTTTGCCTCTGGGATCATCATGGACTCCAGCAT
CTCGAAGCAGGCCCTCAGTGAGATCGAGACCAGGCACAGTGAGATCATCAAGTTGGAGAACAGCATCCGG
GAGCTACAGGATATGTTTCATGGACATGGCCATGCTGGTGGAGAGCCAGGGGGAGATGATTGACAGGATCG
AGTACAATGTGGAACACGCTGTGGACTACGTGGAGAGGGCCGTGTCTGACACCAAGAAGGCCGTCAGTA
CCAGAGCAAGGCACGCAGGAAGAAGATCATGATCATCATTGCTGTGATTCTGGGCATCATCATCGCC
TCCACCATCGGGGCATCTTTGGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >RR206837 representing NM_053788
 Red=Cloning site Green=Tags(s)

MKDRTQELRTAKDSDDDDDVTVTVD RDRFMDEFFEQVEEIRGFIDKIAENVEEVKRKHSAILASPNPDEK
 TKEELEELMSDIKKTANKVRSKLSIEQSIQE EELNRSSADLRIRKTHSTLSRKFVEVMSEYNATQSD
 YRERCKGRIQRQLEITGRITTTSEELEDMLESGNPAIFASGIIMDS SISKQALSEIETRHSEIIKLENSIR
 ELHDMFMMDMAMLVESQGEMIDRIEYNVEHAVDYVERAVSDTKKAVKYQSKARRKIMIIICCVILGIIIA
 STIGGIFG


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_053788

| | |
|-------------------------------|---|
| ORF Size: | 864 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_053788.2 , NP_446240.2 |
| RefSeq Size: | 2130 bp |
| RefSeq ORF: | 867 bp |
| Locus ID: | 116470 |
| UniProt ID: | P32851 |
| Cytogenetics: | 12q12 |
| MW: | 33.1 kDa |
| Gene Summary: | plays a role in neurotransmitter release [RGD, Feb 2006] |