

Product datasheet for **RR206843L4V**

Cxcr5 (NM_053303) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Cxcr5 (NM_053303) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cxcr5
Synonyms:	Blr1; NLR
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_053303
ORF Size:	1122 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206843).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_053303.1 , NP_445755.1
RefSeq Size:	1165 bp
RefSeq ORF:	1125 bp
Locus ID:	29363
UniProt ID:	P34997
Cytogenetics:	8q22
Gene Summary:	G-protein coupled receptor that may play a role in the nervous and immune systems [RGD, Feb 2006]



[View online »](#)