

## Product datasheet for **RR206864L3V**

### Xylt2 (NM\_022296) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Xylt2 (NM_022296) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Xylt2
Synonyms:	xylt-II
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_022296
ORF Size:	2592 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206864).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_022296.2</a> , <a href="#">NP_071632.2</a>
RefSeq Size:	3375 bp
RefSeq ORF:	2595 bp
Locus ID:	64134
Cytogenetics:	10q26
Gene Summary:	isoform of the xylosyltransferase enzyme that catalyzes the transfer of xylose to proteoglycan core protein in glycosaminoglycan biosynthesis [RGD, Feb 2006]



[View online »](#)