

## Product datasheet for **RR206882L3V**

### Commd5 (NM\_139108) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Commd5 (NM_139108) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Commd5
Synonyms:	Hcarg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_139108
ORF Size:	672 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206882).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_139108.2</a> , <a href="#">NP_620808.1</a>
RefSeq Size:	1095 bp
RefSeq ORF:	675 bp
Locus ID:	245974
UniProt ID:	<a href="#">Q9ERR2</a>
Cytogenetics:	7q34
Gene Summary:	plays a role as a regulator of renal epithelial cell growth and differentiation causing G(2)M cell cycle arrest [RGD, Feb 2006]



[View online »](#)