

## Product datasheet for **RR206942L3V**

### Semg1 (NM\_012710) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Semg1 (NM_012710) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Semg1
Synonyms:	RATSVPIIA; Svp1; SVPIIA; Svs2; SVS2P; Svs2p2
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_012710
ORF Size:	1242 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206942).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_012710.2</a> , <a href="#">NP_036842.2</a>
RefSeq Size:	1545 bp
RefSeq ORF:	1245 bp
Locus ID:	25013
Cytogenetics:	3q42
Gene Summary:	comprised of a 13-residue repeat which is involved in the formation of the rat copulatory plug [RGD, Mar 2009]



[View online »](#)