

## Product datasheet for **RR206946L3V**

### Sftpa1 (NM\_017329) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Sftpa1 (NM_017329) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sftpa1
Synonyms:	PSAP; PSP-A; Sftp1; Sftpa; Sftpl; SP-A
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017329
ORF Size:	771 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206946).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017329.2</a> , <a href="#">NP_059025.2</a>
RefSeq Size:	1632 bp
RefSeq ORF:	774 bp
Locus ID:	24773
UniProt ID:	<a href="#">P08427</a>
Cytogenetics:	16p14
Gene Summary:	binds LPS and modifies LPS-induced immune cell activation in the lung [RGD, Feb 2006]



[View online »](#)