

Product datasheet for **RR206948L4V**

Scg2 (NM_022669) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Scg2 (NM_022669) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Scg2
Synonyms:	Chcg
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_022669
ORF Size:	1857 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206948).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_022669.1 , NP_073160.1
RefSeq Size:	2289 bp
RefSeq ORF:	1860 bp
Locus ID:	24765
UniProt ID:	P10362
Cytogenetics:	9q34
Gene Summary:	member of the chromogranin family of proteins that are localized to large dense core vesicles and are endoproteolytically processed to smaller neuropeptides [RGD, Feb 2006]



[View online »](#)