

Product datasheet for **RR206970L4V**

Pou3f3 (NM_138837) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pou3f3 (NM_138837) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pou3f3
Synonyms:	Brn1; RHS1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_138837
ORF Size:	1491 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR206970).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_138837.1 , NP_620192.1
RefSeq Size:	3050 bp
RefSeq ORF:	1494 bp
Locus ID:	192109
UniProt ID:	Q63262
Cytogenetics:	9q22
Gene Summary:	acts as a transcriptional activator in oligodendrocytes; may play a role in regulation of neural development [RGD, Feb 2006]



[View online »](#)