

Product datasheet for RR206982L4V

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

Cd63 (NM_017125) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type: Lentiviral Particles

Product Name: Cd63 (NM_017125) Rat Tagged ORF Clone Lentiviral Particle

Symbol: Cd63

Synonyms: MGC72893

Mammalian Cell Puromycin

Selection: Vector:

pLenti-C-mGFP-P2A-Puro (PS100093)

Tag: mGFP

ACCN: NM_017125

ORF Size: 714 bp

ORF Nucleotide

The ORF insert of this clone is exactly the same as(RR206982).

OTI Disclaimer:

Sequence:

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing

variants is recommended prior to use. More info

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression

varies depending on the nature of the gene.

RefSeg: NM 017125.3, NP 058821.1

 RefSeq Size:
 860 bp

 RefSeq ORF:
 717 bp

 Locus ID:
 29186

 UniProt ID:
 P28648

 Cytogenetics:
 7q11

Gene Summary: human homolog facilitates endocytosis of H,K-ATPase beta-subunit; may play a role in protein

trafficking [RGD, Feb 2006]

