

Product datasheet for **RR207008**

Ldhal6b (NM_183334) Rat Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Ldhal6b (NM_183334) Rat Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ldhal6b
Synonyms:	Ldhl
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>RR207008 representing NM_183334 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**

ATGGCCTGGGCTGTAAGCATCCTGAGGGCCAGAAGAAGAAATGGGTACAGTCAGGGTAGATTTCAAGGACT
TGGCGTAGCTTTATGGCCAGTCAGGTGCTGCGTGTCCCTCCCAACCGGGCTGGTGGAGCCTCACCCC
CAGTCCACCATGGCGACTGTGAAGGTAACCTTCTCCAGCACGTCAGCTCCGAAGAACCGTTCTCCAC
AACCAAGTCTCCATTATAGAACTGGATCCGTCGGCATGGCCTGTGCCATCAGCATCATAGCGAAAGGCT
TGACTGATGAGCTTGCTCTTGTTGATAGTAATGAAGGAAAAATGAAGGGTGAGACAATGGATCTCCAGCA
CGGCAGTGTCTTCATGAAAAATGCCAAATATTGTTTGCAGCAAGGATTTCCATGTCACTGCCAACTCTGAG
GTAGTGGTTATCACAGCGGGAGCACGACAAGAGAAAGATGAGACGCGTCTAAATTTAGTCCAGCGAAACG
TATCTATCTTTAAGATAATGTTTGCCAATATTGTCAAGTACAGCCCTCGCTGCAAAATTCATTGTTTC
CAACCCAGTGGACATCCTAACTACGTAGCGTGGAAATGAGTGGGTTTCCAAAAACCGTGTATTGGA
AGTGGCTGTAACCTAGATACAGCTCGTTTCCGTTACATGCTTGGGCAGAGGCTTGGGATCCACTGAAA
GCTGTGATGGGTGGTCTTGGAGAGCATGGAGACTCAAGTGTCCCTGTGTGGAGTGTAAACATCGC
TGGTGTGCCTCTGAAAGAGCTGAATTCAGCTATAGGAACTAACAAAGACCGGAGCAGTGGGGCAATGTC
CACAAAGAAGTCATTGCCAGTGCCTATAAAAATTACAAAATGAAAGGCTATACTTCTGGCCATTGGCC
TGTCTGTAGTGATATTGCAGAGAGTATTCTGAAGAATCTCAGGAAAACGCACCCAGTTTCCACCAAAAT
TAAAGGTTTATACGGAATAAAAGAAGAGGTGTTCTCAGTGTGCCTTGTGTCCTGGGAGAGAATGGCATA
TCTGATATTATAAAGGTGAAGCTGAGCCCTACAGAGGAGGCCAGATGGTAAAGAGCGCAGACACTATT
GGAAAAATCCAGAAAGATATCAAGTTT

AG**GCGACCG**ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >RR207008 representing NM_183334
Red=Cloning site Green=Tags(s)

MAVAVSILRARRRMGTVRVDFKDLGVALWPSQVLRVPPNRAWWSLTPTSTMATVKGKLLQHVSEEPVLH
 NKVSIIGTGSVGMACAISIIAKGLTDELALVDSNEGKMKGETMDLQHGVSFVMKMPNIVCSKDFHVTANSE
 VVVITAGARQEKDETRLNLVQRNVSIFKIMFANIVKYSRCKIIIVSNPVDILTYVAVKLSGFPKNRVIG
 SGCNLDTARFRYMLGQRLGIHSESCHGWVLGEHGDSSVPVWSGVNIAGVPLKELNSAIGTNKDPEQWGNV
 HKEVIASAYKIIKMKGYTSWAIGLSVADIAESILKNLRKTHPVSTKIKGLYGIKEEVFLSVPCLVGENGI
 SDIIKVKLSPTEEAQMVKSADTIWKIQDKIF


SGPTRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_183334

ORF Size:	1146 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_183334.1 , NP_899163.1
RefSeq Size:	1442 bp
RefSeq ORF:	1149 bp
Locus ID:	369018
Cytogenetics:	1q11
MW:	42 kDa