

Product datasheet for **RR207014L1V**

Pdgfra (NM_012802) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Pdgfra (NM_012802) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Pdgfra
Synonyms:	APDGFR; PDGFACE
Mammalian Cell Selection:	None
Vector:	pLenti-C-Myc-DDK (PS100064)
Tag:	Myc-DDK
ACCN:	NM_012802
ORF Size:	3264 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207014).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_012802.1 , NP_036934.1
RefSeq Size:	6532 bp
RefSeq ORF:	3267 bp
Locus ID:	25267
UniProt ID:	P20786
Cytogenetics:	14p11
Gene Summary:	acts as a receptor tyrosine kinase for PDGF; may play a role in glial cell generation [RGD, Feb 2006]



[View online »](#)