

Product datasheet for **RR207048L4V**

Sdhaf2 (NM_001008371) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Sdhaf2 (NM_001008371) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Sdhaf2
Synonyms:	RGD1309216
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008371
ORF Size:	492 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207048).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001008371.1 , NP_001008372.1
RefSeq Size:	1890 bp
RefSeq ORF:	495 bp
Locus ID:	361726
UniProt ID:	Q5RJQ7
Cytogenetics:	1q43



[View online »](#)

Gene Summary:

Plays an essential role in the assembly of succinate dehydrogenase (SDH), an enzyme complex (also referred to as respiratory complex II) that is a component of both the tricarboxylic acid (TCA) cycle and the mitochondrial electron transport chain, and which couples the oxidation of succinate to fumarate with the reduction of ubiquinone (coenzyme Q) to ubiquinol. Required for flavinylation (covalent attachment of FAD) of the flavoprotein subunit SDHA of the SDH catalytic dimer.[UniProtKB/Swiss-Prot Function]