

## Product datasheet for **RR207076**

### Htr5a (NM\_013148) Rat Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Htr5a (NM\_013148) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Htr5a  
**Synonyms:** 5HT5A; MR22  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR207076 representing NM\_013148  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGATCTGCCTATAAACTTGACCTCCTTTTCTCTCTACTCCCTCCACTTTGGAACCTAACCGCAGCT  
TGGACACGGAAGCCCTGCGCACTAGTCAGTCTTTTCTCTCAGCTTCCGAGTGCTAGTCTGACTTTGCT  
GGGCTTTCTAGCTGCCGCCAGTTCACCTTGAACCTGCTGGTGTGGCCACCATCCTCAGGGTACGCACC  
TTCCACCGAGTACCACACAACCTGGTAGCATCCATGGCTATCTCGGATGTGCTAGTAGCTGTGCTGGTTA  
TGCCCTCTGAGCCTGGTACATGAAGTGTCTGGGCGCCGCTGGCAGCTGGGCGGCGTCTATGCCAGCTGTG  
GATTGCGTGTGACGTCCTCTGCTGTACTGCCAGCATCTGGAATGTGACAGCAATAGCTTTGGACCGCTAC  
TGGTCAATAACGCGCCACCTGGAGTACACACTCCGTGCCCGCAAGCGTGTCTCCAACGTGATGATCCTGC  
TCACTTGGGCACTCTCCGCTGTCATCTCTGGCTCCGCTGCTCTTTGGCTGGGAGAGACTTACTCGGA  
GCTCAGTGAAGAATGCCAGGTGAGTCGCGAGCCTTCTACACCGTGTCTCCACTGTGGGCGCCTTCTAC  
CTGCCGCTGTGTGGTGTCTTTGTATACTGGAAGATTTACAAGGCTGCGAAGTCCGCATGGGCTCCA  
GGAAGACCAACAGCGTCTCCCCATACCTGAAGCTGTGGAGGTGAAGGACGCTTCAACATCCCCAGAT  
GGTGTTCAGTGTCCGTCACGCCACCGTACCTTCCAGACAGAAGGGGACACGTGGAGGGAACAGAAGGAG  
CAAAGGGCAGCCCTCATGGTGGGCATCCTCATTGGGGTGTGCTCTGCTGGTTCCTTCTTTGTTA  
CGGAGCTCATCAGTCCGCTGTGCTCCTGGGATATCCCTGCCCTCTGGAAGAGCATCTTCTGTGGCTGGG  
CTATTCAACTCCTTCTCAACCCGCTCATCTACACGGGTTCAACAGGAGCTACAGCAGTGCTTTTAAG  
GTCTTCTCTCCAAGCAGCAG

**ACGCGT**ACGCGGCGGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR207076 representing NM\_013148  
 Red=Cloning site Green=Tags(s)


MDLPINLTSFSLSTPSTLEPNRSLDTEALRTSQSFLSAFRVLVLTLLGFLAAATFTWNLLVLATILRVRT  
 FHRVPHNLVASMAISDVLVAVLVMPLSLVHELSGRRWQLGRRLCQLWIACDVLCTASIWNVTIALDRY  
 WSITRHLEYTLRARKRVSVMILLTVALSAVISLAPLLFGWGETYSELSEECQVSREPSYTVFSTVGAFY  
 LPLCVVLFVYWKIYKAAKFRMGSRKTNVSPPIPEAVEVKDASQHPQMVFTVRHATVTFQTEGDTWREQKE  
 QRAALMVGILIGVFVLCWFPFFVTELSPLCSWDIPALWKSIFLWLGYNSFFNPLIYTFNRSYSYSAFK  
 VFFSKQQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_013148

<b>ORF Size:</b>	1071 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_013148.1</a> , <a href="#">NP_037280.1</a>
<b>RefSeq Size:</b>	1954 bp
<b>RefSeq ORF:</b>	1074 bp
<b>Locus ID:</b>	25689
<b>UniProt ID:</b>	<a href="#">P35364</a>
<b>Cytogenetics:</b>	4q11
<b>MW:</b>	40.7 kDa
<b>Gene Summary:</b>	receptor for serotonin; involved in mediating serotonergic activity in the brain [RGD, Feb 2006]