

Product datasheet for **RR207133**

Bcl2a1 (NM_133416) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Bcl2a1 (NM_133416) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Bcl2a1
Synonyms: Bcl2a1d
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207133 representing NM_133416
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACAGACTGTGAGTTCATGTATATCCACTCCCTGGCTGAGAACTATCTTCAGTATGTCCTGCAGGTAC
CTGCCTTTGAATCGGCTCCAAGCAAAACGTCCAGAGTGCTACAGAGAGTTGCTTTCTCTGTACAAAAGGA
AGTTGAAAAGAATCTGAAGCCATACTTGGATGACTTTCACGTGGAATCCATAGATACTGCCAGAATAATA
TTCAACCAAGTGATGAAAAAGAATTTGAAGATGGCATCATAACTGGGAAGGATTGTGACTATATTTG
CCTTTGGGGTGTCTCCTGAAAAAGCTTCCACAAGAGCAGATTGGCCTGGATGTGGATACTTACAAGCA
AGTTTCCAGTTTTGTGGCGGAATTCATAATGAATAACACAGGAGAATGGATACAGCAGAATGGAGGCTGG
GAAGATGGCTTCAAAAGAAGTTTGAACCTAAATCTGGCTGGCTGACTTTTCTGCAGATGACAGGGAAGA
TCTGGAAATGCTCTTTCTCCTCAAGCAACTACT

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR207133 representing NM_133416
Red=Cloning site Green=Tags(s)

MTDCEFMYIHLAENYLQYVLQVPAPFESAPSKTSRVLQRVAFSVQKEVEKNLKPYLDDFHVESIDTARII
FNQYMEKEFEDGIINWGRIVTIFAFGGVLLKLPQEQIGLDVDTYKQVSSFVAEFIMNNTGEWIQQNGW
EDGFTKKFEPKSGWLTFLQMTGKIWEMLFLLKQHY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites: SgfI-MluI



Cloning Scheme:



Plasmid Map:



ACCN: NM_133416

ORF Size: 525 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_133416.1</u> , <u>NP_596907.1</u>
RefSeq Size:	699 bp
RefSeq ORF:	528 bp
Locus ID:	170929
Cytogenetics:	8q31
MW:	20.4 kDa