

Product datasheet for **RR207143L4V**

Gpc1 (NM_030828) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Gpc1 (NM_030828) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Gpc1
Synonyms:	HSPG M12; HSPGM12
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_030828
ORF Size:	1674 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207143).
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_030828.1 , NP_110455.1
RefSeq Size:	3509 bp
RefSeq ORF:	1677 bp



[View online »](#)

Locus ID: 58920

UniProt ID: [P35053](#)

Cytogenetics: 9q36

Gene Summary: glycosylphosphatidylinositol-linked heparan sulfate proteoglycan that may play a role in cell adhesion and growth factor response in the nervous system [RGD, Feb 2006]