

## Product datasheet for **RR207147**

### **C1qtnf1 (NM\_001007675) Rat Tagged ORF Clone**

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** C1qtnf1 (NM\_001007675) Rat Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** C1qtnf1  
**Synonyms:** Ctrp1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >RR207147 representing NM\_001007675  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGGCTCCCGTGGACAGGGATTCATGCTGGGATGCTGCCTGCTGCTGGCTGTCACCTGGGCCCCGGTCC  
TGAGCCTTGTGCCACGAGTTCAGGAGGAACAACAGGAGTGGGAAGAGACAGAGGAGCTGACATCTCCCT  
GGATCATGTGACAAGGCCTGAAGAAACACACAAGAGGTATAGCCCTAGCCTGCGTGAGGGCCTCCCCT  
TCCCGGTGCTATCGATGCTGTGATTCCAACCCTCCAGTATACCAGACAATCCCTCCGCCAGATCAACA  
TCACCATCCTAAAAGGTGAGAAAGGAGACCGAGGGGATCGAGGCCTGCAGGGGAAGTATGGCAAAATAGG  
TTCTACGGGTCCCAGGGCCACATTGGCCCCAAAGGGCAGAAGGGATCTATGGGAGCCCCTGGGGACCAC  
TGCAAGAGCCAGTACGCAGCCTTCTGTGGGCCGAAAGAAGGCTCTGCACAGTAACGACTACTTCCAGC  
CTGTGGTCTTCGACACGGAGTTTGTGAACCTCTACAATCACTTCAACATGTTACCCGGAAGTTCTACTG  
CTACGTGCCCGGCATCTACTTCTCAGCCTCAACGTGCACACGTGGAACCAGAAGGAGACCTACCTGCAC  
ATCATGAAGAACGAGGAGGAGGTCGTGATCCTGTATGCGCAGGTGAGCGACCGCAGCATCATGCAGATC  
AGAGCCTGATGCTGGAGCTGCGAGAGCAGGATGAGGTATGGGTGCGCCTTTCAAGGGCGAGCGTGAGAA  
CGCCATCTTCAGTGACGAGTTCGACACCTACATCACCTTCAGCGGCTACCTGGTCAAGCCGGCCTCCGAG  
CCC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >RR207147 representing NM\_001007675  
Red=Cloning site Green=Tags(s)

MGRSGQGFM LGCLLLAVTWGPVLSLVPRVQEEQQEWEETEELT SPLDHSVTRPEETHKRYSPSLREGLPT  
 SRCYRCCDSNPPVYQTI PPPQINITILKGEKGDGRGLQGYGKIGSTGPRGHIGPKGQKGSMPGAPGDH  
 CKSQYAAF SVGRKKALHSNDYFQPVVFDTEFVNLYNHFNMTGKFYCYVPGIYFFSLNVHTWNQKETYLH  
 IMKNEEEVVILYAQVSDRSIMQSQSLMLELREQDEVVWRLFKGERENAI F SDEFDTYIT FSGYL VKPASE  
 P

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001007675

**ORF Size:** 843 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001007675.1</a> , <a href="#">NP_001007676.1</a>
<b>RefSeq Size:</b>	1561 bp
<b>RefSeq ORF:</b>	846 bp
<b>Locus ID:</b>	303701
<b>UniProt ID:</b>	<a href="#">Q5XIG2</a>
<b>Cytogenetics:</b>	10q32.3
<b>MW:</b>	32 kDa