

Product datasheet for **RR207156L4V**

Zhx1 (NM_133620) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Zhx1 (NM_133620) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Zhx1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_133620
ORF Size:	2619 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207156).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_133620.1 , NP_598304.1
RefSeq Size:	4539 bp
RefSeq ORF:	2622 bp
Locus ID:	171159
UniProt ID:	Q8R515
Cytogenetics:	7q33
Gene Summary:	may bind the activation domain of the A subunit of nuclear factor-Y to regulate transcription; localized to the nucleus; can form a homodimer in vitro [RGD, Feb 2006]



[View online »](#)