

## Product datasheet for **RR207163L3V**

### Ntm (NM\_017354) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Ntm (NM_017354) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Ntm
Synonyms:	Hnt; RNU16845
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_017354
ORF Size:	1032 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207163).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_017354.1</a> , <a href="#">NP_059050.1</a>
RefSeq Size:	2040 bp
RefSeq ORF:	1035 bp
Locus ID:	50864
UniProt ID:	<a href="#">Q62718</a>
Cytogenetics:	8q13
Gene Summary:	member of a family of GPI-anchored cell adhesion molecules; may be involved in specifying neuronal projections and connectivity [RGD, Feb 2006]



[View online »](#)