

## Product datasheet for **RR207192L3V**

### **Rps6ka5 (NM\_001108048) Rat Tagged ORF Clone Lentiviral Particle**

#### **Product data:**

Product Type:	Lentiviral Particles
Product Name:	Rps6ka5 (NM_001108048) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Rps6ka5
Synonyms:	Msk1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001108048
ORF Size:	1554 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207192).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001108048.1</a> , <a href="#">NP_001101518.1</a>
RefSeq Size:	5427 bp
RefSeq ORF:	1557 bp
Locus ID:	314384
Cytogenetics:	6q32



[View online »](#)