

Product datasheet for **RR207229L4V**

PPP3R1 (NM_017309) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	PPP3R1 (NM_017309) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	PPP3R1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017309
ORF Size:	510 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207229).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017309.2 , NP_059005.1
RefSeq Size:	2803 bp
RefSeq ORF:	513 bp
Locus ID:	29748
UniProt ID:	P63100
Cytogenetics:	14q22
Gene Summary:	regulatory subunit of the enzyme that catalyzes the calcium-dependent removal of serine- and threonine-bound phosphate groups [RGD, Feb 2006]



[View online »](#)