

Product datasheet for **RR207230**

Cd27 (NM_001024335) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Cd27 (NM_001024335) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Cd27
Synonyms: Tnfrsf7
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207230 representing NM_001024335
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCATGGCCACCCCTCTACTGGCTCTGCATGCTGGGACCTTGGTAGGGCTCTTAGCTACCCAGCCC
CAAACAACGTCCAGACAGACACTACTGGATTGGGGCAGGACTCTGCTGCCAGATGTGTGGCCAGGGAC
ATTCCTTGTGAAGCACTGTGACCAGGACAGAGCCGCTGCTCAGTGTGATCCGTGTATTCCAGGCAGTCC
TTCTCTCCAGACTACCACACCCGGCCCACTGCGAGAGCTGCAGGCACTGTAACCTCTGGTTTTCTCATCC
GCAACTGCACGGTCACTGCCAATGCTGAGTGCACCTGTTCCAAGGGCTGGCAGTGCAGGGACCAGGAATG
TACGGAGTGTGACCCTCCTCTAAACCCTGCACTGACTAGCCAGCCCTCTGAGGCCCCAGCCCCAGCTG
CCACCACCCACCACTTACCTTATGCCACAGAGAAGCCATCCTGGCCCCACAAAGGCAGCTTCCCGACT
CAACTGTCTATAGCCGCTGCCATCCAAAGACCTCTGTGCAGCTCCGACTGTATCCGGATCTTTGTGAC
CTTCTCCAGCATGCTTCTGTCTTCGTCCTGGGTGGGATCTTGTCTTCCATCAAAGAAGAAATCATGGG
CCAAATGAAGATAGCCAGGCAGTACCTGAAGAGCTTTGTCCTTACAGCTGCCCTAGGGAGGAGGAGGGCA
GTGTCATCCCTATCCAGGAAGACTACCGAAACCTGAGCCTGCTTCTACCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >RR207230 representing NM_001024335
 Red=Cloning site Green=Tags(s)

MAWPPLYWLCMLGTLVGLLATPAPNNPCDRHYWIGAGLCCQMCGPGTFLVKHCDQDRAAAQDCPCIPGTS
 FSPDYHTRPHCESCRHCNSGFLIRNCTVTANAECTCSKGWQCRDQECTEDPPLNPALTSQPSEAPSPQL
 PPPHTLPYATEKPSWPPQRQLPDSVYSLRPSQRPLCSSDCIRIFVTFSSMLLVFLVGGILFFHQRRNHG
 PNEDSQAVPEELCPYSCPREEEGVIPIQEDYRKPEPASYP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001024335
ORF Size: 753 bp

OTI Disclaimer: Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001024335.1](#), [NP_001019506.1](#)

RefSeq Size: 1602 bp

RefSeq ORF: 756 bp

Locus ID: 500318

Cytogenetics: 4q42

MW: 27.9 kDa