

Product datasheet for **RR207249**

Reep5 (NM_001108888) Rat Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Reep5 (NM_001108888) Rat Tagged ORF Clone
Tag: Myc-DDK
Symbol: Reep5
Synonyms: Dp1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >RR207249 representing NM_001108888
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCCGACCCATGAGGGAGAGGTTTCGACCGGTTCTGCACGAGAAGAACTGCATGACTGATCTCCTGG
CCAAGCTCGAGGCCAAGACTGGAGTGAACCGGAGCTTCATCGCGCTCGGTGTCATCGACTGGTGGCTTT
GTATCTGGTGTTCGGTTATGGAGCCTCTCTCTCTGCAACCTGATAGGTTTCGGATACCCAGCCTACATC
TCAATGAAAGCCATCGAGAGCCCAACAAGATGATGACACCCAGTGGCTGACGTACTGGTGGTGTATG
GTGTGTTGAGCATCGCCGAGTCTTCTCCGACCTTTCTCTGTCCTGGTCCCTTCTACTACATGCTGAA
GTGCGGCTTCTGTTGTGGTGCATGGCTCCCAGCCGTCACCGGGCTGAACTGCTGTACCGGCGTGTGTC
ATCCGTCCTTCTCTGAAGCACGAGTCCCAGGTAGACAGTGTGGTGAAGGACGTGAAGGACAAAGCCA
AAGAGACTGCAGACGCCATCAGCAAAGAAGTCAAGAAGGCTACGGTGAACCTGCTGGCGCAGAGAAGAA
GAGCACC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >RR207249 representing NM_001108888
Red=Cloning site Green=Tags(s)

MSAAMRERFDRFLHEKNCMTDLLAKLEAKTGVNRSFIALGVIGLVALYL VFGYGASLLCNLIGGYPAYI
SMKAIESPNKDDDTQWLYWVYGVFSIAEFFSDLFLSWFPFYMLKCGFLLWCMAPSPSNGAELLYRRY
IRPIFLKHESQVDSVVKDVKAKETADAISKEVKKATVNLLGDEKKST

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI




[View online »](#)

Cloning Scheme:



Plasmid Map:



ACCN: NM_001108888

ORF Size: 567 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001108888.2](#), [NP_001102358.2](#)

RefSeq Size: 2917 bp

RefSeq ORF: 570 bp

Locus ID: 364838

UniProt ID: [B2RZ37](#)

Cytogenetics: 18p12

MW: 21.4 kDa

Gene Summary: May promote functional cell surface expression of olfactory receptors.[UniProtKB/Swiss-Prot Function]