

Product datasheet for **RR207249L3V**

Reep5 (NM_001108888) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Reep5 (NM_001108888) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Reep5
Synonyms:	Dp1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-Myc-DDK-P2A-Puro (PS100092)
Tag:	Myc-DDK
ACCN:	NM_001108888
ORF Size:	567 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207249).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_001108888.2 , NP_001102358.2
RefSeq Size:	2917 bp
RefSeq ORF:	570 bp
Locus ID:	364838
UniProt ID:	B2RZ37
Cytogenetics:	18p12
Gene Summary:	May promote functional cell surface expression of olfactory receptors.[UniProtKB/Swiss-Prot Function]



[View online »](#)