

## Product datasheet for **RR207272L4V**

### Cldn19 (NM\_001008514) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Cldn19 (NM_001008514) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Cldn19
Synonyms:	claudin-19
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_001008514
ORF Size:	633 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207272).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_001008514.1</a> , <a href="#">NP_001008514.1</a>
RefSeq Size:	1626 bp
RefSeq ORF:	636 bp
Locus ID:	298487
UniProt ID:	<a href="#">Q5QT56</a>
Cytogenetics:	5q36
Gene Summary:	Plays a major role in tight junction-specific obliteration of the intercellular space, through calcium-independent cell-adhesion activity.[UniProtKB/Swiss-Prot Function]



[View online »](#)