

Product datasheet for **RR207318L4V**

Slc7a5 (NM_017353) Rat Tagged ORF Clone Lentiviral Particle

Product data:

Product Type:	Lentiviral Particles
Product Name:	Slc7a5 (NM_017353) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Slc7a5
Synonyms:	E16; TA1
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_017353
ORF Size:	1536 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207318).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	NM_017353.1 , NP_059049.1
RefSeq Size:	3455 bp
RefSeq ORF:	1539 bp
Locus ID:	50719
UniProt ID:	Q63016
Cytogenetics:	19q12
Gene Summary:	integral membrane proteininvolved in liver development, carcinogenesis, and cell activation [RGD, Feb 2006]



[View online »](#)