

## Product datasheet for **RR207323L4V**

### Dpysl5 (NM\_023023) Rat Tagged ORF Clone Lentiviral Particle

#### Product data:

Product Type:	Lentiviral Particles
Product Name:	Dpysl5 (NM_023023) Rat Tagged ORF Clone Lentiviral Particle
Symbol:	Dpysl5
Synonyms:	Crmp5; Ulip6
Mammalian Cell Selection:	Puromycin
Vector:	pLenti-C-mGFP-P2A-Puro (PS100093)
Tag:	mGFP
ACCN:	NM_023023
ORF Size:	1692 bp
ORF Nucleotide Sequence:	The ORF insert of this clone is exactly the same as(RR207323).
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
RefSeq:	<a href="#">NM_023023.1</a> , <a href="#">NP_075412.1</a>
RefSeq Size:	4845 bp
RefSeq ORF:	1695 bp
Locus ID:	65208
UniProt ID:	<a href="#">Q9JHU0</a>
Cytogenetics:	6q14
Gene Summary:	may have a role in neuronal differentiation and/or axon growth [RGD, Feb 2006]



[View online »](#)